


Application list No.	38
Title	Determination of bromates in water samples
Types of samples	All kinds of disinfected water samples
Principle	<p>The bromate anions in the stabilised water sample are reduced to bromine with an acidified bromide solution:</p> $\text{BrO}_3^- + 5 \text{Br}^- + 6 \text{H}^+ \rightleftharpoons 3 \text{Br}_2 + 3 \text{H}_2\text{O}$ <p>The formed bromine is then determined by in-electrode coulometric titration through its electrochemical reduction to bromide ions:</p> $\text{Br}_2 \rightleftharpoons 2 \text{Br}^- - 2 e^-$ <p>During this step the chronopotentiometric signal is recorded and evaluated. Possible interfering species are removed by boiling and addition of Fe(II) ions.</p>
Reagents	<p>§ HCl, p.a., conc.</p> <p>§ FeSO₄·7 H₂O, p.a.</p> <p>§ Certified reference material for bromate, e.g. dried KBrO₃ or NaBrO₃</p>
Solutions	<p>§ R-013 Carrier electrolyte, 0,1 mol/dm³ HCl</p> <p>§ R-025 Ethylenediamine (EDA) in water, 100 mg/cm³</p> <p>§ R-026 (Fe²⁺ solution: 0.5 g FeSO₄·7 H₂O dissolved in 100 ml 0.3 mol/dm³ H₂SO₄)</p> <p>§ R-027 (KBr solution)</p> <p>§ Bromate standard solution in water (50 mg/dm³ BrO₃⁻)</p>
Standard solutions for the calibration curve	<p>The calibration solutions are prepared from the bromate standard solution (see above) in water containing 50 mg/dm³ EDA::</p> <p>Standard No. 1: 5 µg/dm³ BrO₃⁻</p> <p>Standard No. 2: 10 µg/dm³ BrO₃⁻</p> <p>Standard No. 3: 20 µg/dm³ BrO₃⁻</p> <p><i>Note: If appropriate, adjust the composition of the standard solutions to the sample!</i></p> <p>The shelf life of the calibration solutions is one day.</p> <p>The calibration solutions are treated in the same way as the samples – see the sample preparation procedure.</p>
Sampling and sample preservation	<p><i>U.S. EPA Method 300.1:</i></p> <ol style="list-style-type: none"> The water samples are sampled into amber glass flasks. Add ethylenediamine (R-025) so that the final concentration is 50 mg/dm³. The samples can be stored at or below 6 °C but not longer than 28 days.
Sample preparation	<ol style="list-style-type: none"> Boil samples containing disinfection byproducts such as chlorine, hypochloride, chloramines, ozone, chlorine dioxide, chlorates, as well as carbon dioxide for 5 min and let it cool to laboratory temperature. To a dry 100 ml volumetric flask add 0.5 ml of Fe(II) solution (R-026) (1.0 ml for samples with HCO₃²⁻ contents over 500 mg/l) Adjust the volume to 100 ml with the sample (or with the calibration solution in the case of a calibration) and on mixing let the solution stay for 2 to 10 min. Important: the pH must lay in the range of 5 - 6! (In the case of the calibration solutions it may be pH 3 to 6) Add 2 ml of concentrated HCl and mix the solution: Solution A To a dry 50 ml volumetric flask add 0.1 ml of the R-027 reagent and adjust the volume to 50 ml with the solution A. Mix the solution: Solution B. Solution A serves as the blank for the treated sample (or standard). Solution B serves as the sample (or standard). These solutions are stable not longer than 30 min and therefore must be immediately analysed.
Electrode	E104 C (moist the new electrode with ethanol before setting into the cell)
Experimental parameters	Galvanostatic mode

	<div style="border: 1px solid black; padding: 5px;"> <div style="display: flex; justify-content: space-between;"> <div style="width: 30%;"> <p>Calibration</p> <input type="radio"/> Calibrationless <input checked="" type="radio"/> Calibration curve <input type="radio"/> Standard addition </div> <div style="width: 30%;"> <p>Deposition</p> <input checked="" type="radio"/> GST <input type="radio"/> PST </div> <div style="width: 30%;"> <p>Background reading</p> <input type="radio"/> With the first measurement only <input checked="" type="radio"/> With each measurement <input type="radio"/> With each new sample or standard <input type="radio"/> With each <input type="text" value="0"/> th sample or standard </div> </div> <div style="border: 1px solid black; padding: 5px; margin-top: 5px;"> <p style="text-align: center; background-color: #0056b3; color: white; margin: 0;">Setup Parameters</p> <div style="display: flex; justify-content: space-between; border-bottom: 1px solid black; margin-bottom: 5px;"> General Preparation Regeneration Measure Calibration Calculation Samples </div> <div style="display: flex; justify-content: space-between;"> <div style="width: 30%;"> <p>Deposition</p> <p>Edepos: <input type="text" value="600"/> mV</p> <p>Idepos: <input type="text" value="0"/> uA</p> </div> <div style="width: 30%;"> <p>Pause (s)</p> <p>Quiesc1: <input type="text" value="0"/></p> <p>Quiesc2: <input type="text" value="0"/></p> <p>Regen: <input type="text" value="1"/></p> </div> <div style="width: 30%;"> <p>Flow (ml)</p> <p>Sample: <input type="text" value="4"/></p> <p>Back: <input type="text" value="4"/></p> <p>Rinse: <input type="text" value="0"/></p> </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 30%;"> <p>Potencial (mV)</p> <p>Estart1: <input type="text" value="550"/></p> <p>Estart2: <input type="text" value="550"/></p> <p>Estop: <input type="text" value="100"/></p> <p>Eregen: <input type="text" value="300"/></p> <p>Estdby: <input type="text" value="300"/></p> </div> <div style="width: 30%;"> <p>Stripping:</p> <p>Istrip: <input type="text" value="-10"/> uA</p> <p>Timeout: <input type="text" value="180"/> s</p> <p>Pump: <input checked="" type="radio"/> Off <input type="radio"/> On</p> </div> <div style="width: 30%;"> <p>Sample segmentation</p> <input checked="" type="radio"/> No <input type="radio"/> Yes </div> </div> <p style="text-align: right; margin-top: 5px;">Flow: <input type="text" value="6"/> ml/min</p> </div> </div>
<p><i>Typical signal</i></p>	
<p><i>Data evaluation</i></p>	<p>Calibration curve</p>
<p><i>Metrological data</i></p>	<p>Concentration range: 0.5 to 50 µg/dm³ BrO₃⁻</p> <p>Reproducibility: 6,8 % at 5 µg/dm³ BrO₃⁻ in the measured solution 3,1 % at 20 µg/dm³ BrO₃⁻ in the measured solution</p>
<p><i>Interferences</i></p>	<p>§ High Fe(III) contents (over 10 mg/l) deteriorate the results. Use a cation exchanger to remove Fe species.</p> <p>§ High nitrate contents (over 50 mg/dm³) interfere</p>
<p><i>Notes</i></p>	<p>§ Samples with bromate contents over 50 µg/dm³ should be diluted with water prior to analysis.</p> <p>§ The reagent solutions should be stored in a refrigerator.</p> <p>§ The analyses can be accelerated by enhancing the flow rate to 6 ml/min through replacing the pump tube.</p> <p>§ If on addition of the Fe(II) reagent to the sample solution the resulting solution gains a yellowish colour, it indicates a pH value of the resulting solution over 6 and the analysis would probably be wrong. In such a case adjust the sample pH, e.g. by adding some diluted HCl, so that on adding the Fe(II) reagent the resulting solution would have a pH in the range of 5-6: the solution should be virtually colourless.</p> <p>§ Removal of air bubbles from the system:</p> <ul style="list-style-type: none"> - Immerse the sampling tube to ethanol - Start the „Preparation“ routine - Make 2-3 measurement of a standard solution, e.g. 10 µg/l BrO₃⁻.